



Be Right™



Free Ammonia Chemkey Reagents (box of 25)

Product #:	9429500
AED Price:	Contact Hach

Free Ammonia Chemkeys, range 0.05 - 0.5 mg/L NH₃-N. For SL1000 Portable Parallel Analyser (PPA).

Specifications

Instrument:	SL1000 Portable Parallel Analyser
Measuring range:	0.05 - 0.50 mg/L NH ₃ -N
Parameter:	Ammonia, Free
Platform:	Chemkeys
Quantity:	pk/25