



Ammonium cuvette test 0.005-0.05 mg/L $\text{NH}_4\text{-N}$, 20 tests

Product #:

LCK504

AED Price:

Contact Hach

Hazardous



LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Environmentally friendly due to less chemicals.

Convenient handling with 50 mm rectangular cuvettes

Highest accuracy for analysis of ultra low range concentrations

Specifications

According to standard:	ISO 7150-1, DIN 38406 E5-1, UNI 11669:2017, ISO 23695
Instrument:	DR3900, DR6000
Measuring range:	0.005 - 0.05 mg/L $\text{NH}_4\text{-N}$
Measuring range (2):	0.006 - 0.06 mg/L $\text{NH}_4\text{-N}$
Method:	Indophenol blue
Number of tests:	20
Parameter:	Ammonium
Platform:	LCK
Storage conditions:	2 - 8 °C (keep refrigerated)