



Phosphate (ortho/total) cuvette test 10-100 mg/L PO₄-P, 25 tests

Product #:

LCK351

AED Price:

Contact Hach

Hazardous



LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Officially approved for consent limits. Environmentally friendly due to less chemicals.

Precise and reliable measurement values

Maximum safety for users, thanks to the closed system and low amounts of reagents

Convenient and error-free dosing of the reagents

Barcode label for automatic recognition in the photometer and Truecal calibration

Differentiation between tests and measuring ranges by means of colour-coding

Specifications

According to standard:	2321142
Description:	Ortho/Total Phosphate
Digestion required:	for Total P
EPA compliant:	N/A
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Measuring range:	10 - 100 mg/L PO ₄ -P
Measuring range (2):	30 - 300 mg/L PO ₄
Method:	Phosphormolybdenum Blue
Number of tests:	25
Parameter:	Phosphate, ortho + total
Platform:	LCK
Standard method:	EN ISO 6878
Storage conditions:	15 - 25 °C
Truecal:	Yes

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)