



Total Hardness Reagent Set for Digital Titrator (LR)

Product #:

2448000-EU

AED Price:
Available

Contact Hach

Hazardous



For the determination of Total Hardness, 10-160 mg/L CaCO_3 , EDTA Method.

Specifications

EPA compliant:	Yes
Measuring range:	10 - 160 mg/L CaCO_3
Method:	EDTA Titration
Method number:	8213
Number of tests:	40-100
Parameter:	Total Hardness
Platform:	Digital Titrator
Special Features:	Reagent set, LR, Digital Titrator

Required Accessories

- Protective gloves size L, blue, nitrile, powder-free, 50 pieces (Item SM743L)