



NitraVer 5 Nitrate Reagent Powder Pillows, 5 mL, pk/100

Product #:

1403599-EU

Hazardous

AED Price:
Available

Contact Hach



Powder reagent in sealed packages for the determination of Nitrogen-Nitrate up to 50 mg/l NO₃-N.

For high range nitrate determination by the NitraVer 5 cadmium reduction method. Range: up to 50 mg/L as nitrate-N. Sample Size: 5 mL. Pack of 100 powder pillows.

Specifications

Method:	Cadmium Reduction
Platform:	Powder Pillows
Quantity:	100 /pk
Sample size:	5 mL

Required Accessories

- Protective gloves size L, blue, nitrile, powder-free, 50 pieces (Item SM743L)