



Lead reagent set, 5-150 µg/L Pb

Product #: 2375000-EU
AED Price: Contact Hach
Available

Hazardous



Acid soluble Lead in a potable water sample is first concentrated on a Fast Column Extractor. The Lead is then eluted from the Extractor and determined colorimetrically with an indicator.

Specifications

Instrument:	See Methods/Procedures for instruments compatible with this product
Measuring range:	5 - 150 µg/L Pb
Method:	LeadTrak
Method number:	8317
Number of tests:	20
Parameter:	Lead
Shipping Method:	Lead Reagent set

Required Accessories

- Protective gloves size L, blue, nitrile, powder-free, 50 pieces (Item SM743L)